

Copy of HIRES Setup Instructions to be sent to Keck Support Astronomers (8 Feb 2018):

Setting Up HIRES

- Check WEB-based Status of HIRES before starting. When HIRES is released, you may begin.
 - <http://www2.keck.hawaii.edu/realpublic/inst/siasrng/release/web/Sias/siasEh.php?host=mysqlserver&table=datebook>
- From tan background, drag & click on HIRES CONTROL MENU
 - Drag over to: START all HIRES GUIs w/ hexpocon
 - Answer questions:
 - "Do you want to continue running the setup script?" type y and return
 - Enter observer names
 - Confirm data directory (/s/sdata125/hires1/20??monDD/)
 - Use the current UT date. (typically one day later than local calendar date), hit enter.
 - If we have 2nd half ONLY, the append "_B" to the data directory.
 - Set starting observation number at 1 on first night of run.
 - Now appearing on the 3 screens:
 - 1) HIRES dashboard, Exposure meter DashBoard, terminal
 - 2) XHIRES GUI, dewar level window
 - 3) SAO image: ds9
- Check dewar level. In a terminal type: ln2
 - If no terminal is present, right click on tan background and goto: Keck Menu > Xterms > Xterm (needs corrected)
 - The dewar level is also visible when using "START all HIRES GUIs" in the window with the XHIRES GUI.
 - Top off dewar (if level is below 70%) by: right-clicking on tan background,
 - drag to "HIRES Control Menu" and "Initiate HIRES Dewar Fill"
 - dewar evaporation rate is 5% per hour and auto-refills at 10%. Always try avoid an autorefill.
 - If the dewar is filled after 5pm Pacific Time, then it does not need refilled near sunset.
- START THE IODINE CELL
 - Select from desktop pulldown: "HIRES control menu > Iodine cell menu > Start Iodine Cell"
 - Cell takes 45 minutes to warm up fully. A warm cell reads:
 - tempiod1 | 65 degrees
 - tempiod2 | 50 degrees (+/- 0.1 degree)
- Open the Mirror Covers
 - From XHIRES GUI > click ETC > Click OPEN RED
- **First** Night of run: Set up data directories.
 - On HIRES dashboard, click on yellow "Start Here" button.
 - Click on "retrieve" to install directory for raw data and frame number.
 - Update "Filename root" with current run number. (j2??)
 - If second (or later) night, round file number up to next hundred
 - Click on "Commit" to set values.
 - At the command line, create all of the directories needed for the **entire** run.
 - for example: /s/sdata125/hires1/2016nov15/ UT Date
 - If second half only, append "_B" to the directory name for that UT date.
 - Double-click to the right of "FULL FRAME" button on HIRES DashBoard to view CCD control panel.
 - Make sure Binning is set to 3x1 by clicking on binning button : "x3y1"
 - Double check all values in logsheet header are correct.
 - (i.e., copy most recent logsheet file, update to new name (e.g. j239.logsheet1))
- **Second** night of run or later,
 - In HIRES Dashboard, change the "Data directory" by clicking on "User Config"/"Start Here."
 - Update Data Directory with UT date of observation
 - Press 'Commit'
 - Increment the frame number to the next even hundred.
 - Update Observer names
 - Double check all values in logsheet header are correct.

Cut to the Chase Check list

- Once the appropriate GUIs are up and running, step through:
 - All commands must be run in lehoula window.
- check dewar level: ln2
- top off dewar (if necessary):

- click on tan background, drag to HIRES control, drag to "Initiate HIRES Dewar Fill"
- XHIRES GUI:
 - From 'etc' button on xhires upper-left: open-red
 - Slit should read 14.08" (m slitname = opened)
 - Filter1 = clear (m fil1name = clear; using gui okay) (Formerly KV370)
 - Filter2 = clear (m fil2name = clear; using gui okay)
 - Collimator = red (This should always be set by SA)
- In a terminal:
 - Set Collimator focus: m cofraw = +70000 ; Use 's cofraw' to show value
 - Set Camera focus: m cafrw = 0 ; Use 's cafrw' to show value
 - Update Echelle and Cross disperser angles with "A" button
- Next to each element on XHIRES GUI. Use values from a recent logsheet.
 - Set Echelle
 - Set Cross disperser
 - Avoid single movements > 0.5 to help prevent oscillations. Move in multiple steps if needed.
 - If no response, in lehoula: stop_cd_drift
 - If oscillating, from background: HIRES control menu > Stop Cross-disperser Oscillation
- TV filter: BG38 + ND0.01 (BG38 is important, ND up to OA)
- On HIRES Dashboard:
 - Double-Click on brown "CCD" rectangle, right of "FULL FRAME"
 - CCD Binning: Enter in the left box: X = 3 , Y = 1 or click on "Binning" and pull down to "X3Y1".
 - No Need to set any other params in CCD window.
 - Check/Set OUTDIR: directory for raw data
 - Check/Set OUTFILE: prefix of filenames, i.e. j160
 - Check CCD readout mode:
 - Gain = "low" (default)
 - Speed = "fast"
- Check that iodine cell is on and temperature window appears.
 - If trouble, try: setenv DISPLAY lumahai:0.0
 - xshow -s hires tempiod1 tempiod2 &
 - Set ObsType name field on brown CCD rectangle as "Object."

Logsheets

- Copy a recent logsheet to a new one, i.e.,
 - cp j200.logsheet1 j201.logsheet1
 - Update Observers, run number (first night of run)
 - Update UT date (eg 15 May 2011) (NOT 14/15 May 2011)
 - Periodically note the seeing with column cut: FWHM: 0".38/bin
 - This is a check of the seeing measured by eye.

MID-TIME DETERMINATION

- Login to a Keck machine (e.g., on your laptop – instructions above)
- To determine the weighted midpoints, and print an nice line for the logsheet run mtfits.pro
 - cd ~gmarcy/focus (on keck computer)
 - In IDL,
 - IDL> mtfits, inpdire='/s/sdata125/hires1/YYYYmonDD/', 'j420023', n=5
 - Prints logsheet line for j420023.fits and subsequent four exposures.
 - Copy output of MTFITS into logsheet.
 - Put any flux-weighted midpoint error messages in comment field of logsheet entry for that observation.
 - If there are issues running mtfits.pro on the Keck computers, it can also be run from /mir3/focus/ if the raw data is present on cadence.

Focus

Focus using Thorium Lines

- Turn OFF Exposure meter
- Lamp : Th-Ar #2
- Lamp Filter: NG3 filter
- In terminal window:
 - m deckname = D5 (in lehoula window)
- Iodine : Out
- ObsType = Object
- Texp : 10 sec (in brown "CCD" window, enter exposure time. Click "UpdateCCD")
- Click "EXPOSE"
 - **Perform focus in IDL.**
 - cd to focus directory
 - To specify a file for focus:
 - **IDL> foc,inpfile='jnnnxxxx.fits',/mark,/plt**
 - If necessary, compare lines to list included at the end of this document.
 - Check instructions from focus program and move echelle and cross-disperser as needed.
 - If focus program crashes, you may need to move echelle or cross disperser manually.
 - If note regarding 'Counts in lines too low' appears. Re-position lines manually.

Manual Grating/Echelle Move Manually

- Grating Move:
 - Horizontal: +0.001 deg of echelle rot moves lines left by 1 column
 - Vertical: +0.002 deg of X-disp rot moves lines down by 1 row
 - + As a last resort change cafrac or cofrac on command line to focus
 - (using m cafrac= and m cofrac=). Try cafrac first; steps of ~10,000 are needed in cafrac to make any appreciable difference in focus.
 - Record nominal FWHM, echelle position, and X-disperser values in the logsheet.

Calibrations

- THORIUM Exposures w/ B5
 - Turn OFF Exposure meter
 - Lamp : Th-Ar #2
 - Lamp Filter: ng3
 - m deckname = B5 (0.85 x 3.5 arcsec, ==> 4.0 pixel projected slit)
 - WARNING: use m deckname=B5, NOT the HIRES GUI.
 - Iodine : Out
 - Exposure : 1 sec (take 1 or 2 at begin and end of night)
- THORIUM Exposure w/ B1
 - Turn OFF Exposure meter
 - Lamp : Th-Ar #2
 - Lamp Filter: ng3
 - m deckname = B1 (0.57 x 3.5 arcsec, ==> 3.0 pixel projected slit)
 - WARNING: use m deckname=B1, NOT the HIRES GUI.
 - Iodine : Out
 - Exposure : 2 sec (take 1 or 2 at begin and end of night)
- Iodine Cell Calibrations w/ B1
 - Make sure cell is fully warmed up (see p.1) before taking these.
 - Turn OFF Exposure meter.
 - Lamp : Quartz2
 - Lamp Filter: ng3
 - Aperture : B1 (0.85 x 3.5 arcsec, ==> 4.0 pixel projected slit)
 - WARNING: use m deckname=B1, NOT the HIRES GUI.
 - Iodine : In
 - Exposure : 3 secs
 - *check saturation*: < 20,000 counts on middle chip?
 - *Check I2 line depth*. In center of chip, it should be ~30%
 - column cut with DS9: Region > more... > proj
- Iodine Cell Calibrations B5
 - Make sure cell is fully warmed up before taking these.
 - Turn OFF Exposure meter.
 - Lamp : Quartz2

- Lamp Filter: ng3
- Aperture : B5 (0.85 x 3.5 arcsec, ==> 4.0 pixel projected slit)
 - WARNING: use m deckname=B5, NOT the HIRES GUI.
- Iodine : In
- Exposure : 2 secs
 - *Check saturation:* < 20,000 counts on middle chip?
 - *Check I2 line depth.* In center of chip, it should be ~30%
- column cut with DS9: Region > more... > proj

Wide Flat-Fields

- Turn OFF Exposure meter
- Lamp : Quartz2
- Lamp Filter: ng3
- m deckname = C1 (0.85 x 7.0 arcsec, ==> 4.0 pixel projected slit)
 - WARNING: use m deckname=C1, NOT the HIRES GUI.
- Iodine : Out
- Exposure : 1 sec (check saturation: middle chip should have 10,000 < counts < 20,000)
- Take 50 exposures at the beginning of the night
 - Check one test exp for saturation (<20k counts), then take multiple exposures with command:
 - exp_acq 49 ; m lampname=none ; m deckname=C2
 - Concatenates commands to turn off lamp after flats finish.
- Top off LN dewar if necessary.

Observations of Stars

- Iodine : IN (check I2 temp: 50C)
- ~30 min before sunset:
 - top off LN dewar if necessary (definitely refill if level is < 100.0... drops ~5% per hour)
- Just before observing on screen 4(right-most, left-most at CIT):
 - From tan background click and select K1 Guider Eavesdropping > Start Observer UI
 - type y, enter
 - From tan background click and select K1 Telescope Status Menu > FACSUM
 - From tan background click and select K1 Telescope Status Menu > XMET
- Start exposure meter.
 - Click on the upper left button "System Start" on exposure meter.
 - Click on "Arm" in upper left of right panel to start target monitoring.
 - Default exposure level is 250000 equivalent to SNR ~200
- Set max exposure to 500s (in HIRES Dashboard CCD ExpTime)
- Expected Exposure time:: At V=8, S/N=300 in 300 seconds
- Remember to open HIRES hatch
- Check that
- Check with OA that "slit guiding algorithm" is being used.
- Add a "fiducial mark" at the position of the star by right clicking the magic guider snapshot at the desired location
- Double check:
 - that iodine temperature is 50/65C, and iodine is running
 - vertical angle mode is on and set to 0.
 - Filter #1 is "clear"
 - Filter #2 is "clear"
 - TV filters are "bg38" and "nd_0.01"
- Iodine : IN
- In Twilight:

- m deckname = C2 (0.85 x 14.0 arcsec)
- If seeing is > 2.0", then begin observing only 10 minutes before 12deg twilight and use B5.
- Likewise if seeing > 2.0" at the end of the night, use B5 in twilight and end 10 minutes after 12deg
- WARNING: use command line, NOT HIRES GUI
- Non-twilight:
 - m deckname = B5 (0.85 x 3.5 arcsec, ==> 4.0 pixel projected slit)
 - WARNING: use command line, NOT HIRES GUI
- Poor seeing > 2 arcsec
 - Stick to V < 10 stars (throughput)
 - Use B5 decker. Sky subtraction does not work well when stellar PSF fills the slit (seeing > 2.5").
- Telescope wrap limits
 - From the south wrap, moving through the west, the north limit is an azimuth of 325 degrees.
 - From the north wrap, moving through the west, the south limit is an azimuth of 235 degrees.

END of NIGHT Procedure

- Turn off exposure meter.
- Close the hatch
- Take B1/B5 iodine exposures.
- Take B1/B5 thorium exposures.
- Turn off lamps, but leave everything else open, if not last night of run.
- Logsheet:
 - Complete logsheet.
 - Make sure there are no missing or duplicate lines

Shutdown Sequence

- On last night of a run run the following:
 - From background pulldown, HIRES control menu > End of Night Shutdown

Odds and Ends

Logsheet Etiquette

- Change calibration names to 'focus','thar','iodine','wideflat'
- Try to use only these names for calibrations, unless running specific tests.
- Every calibration(except wideflat) needs its own line in the logsheet.

Tips, Tricks, & Troubleshooting

- Cross-disperser Oscillations:
- If cross-disperser values are oscillating, reset by right-clicking blue background and going to HIRES Control Menu > Stop Cross-disperser Oscillation.
- Avoid moving cross-disperser by increments > 0.5 to help prevent oscillations. Move in multiple steps if needed.
- The GUIs.
 - HIRES : DashBoard
 - Controls exposure lengths & file writeouts
- Appears with the "Start all HIRES GUIs" initialization sequence.
- Check dewar levels: ln2.
- Display cofraw, cafraw, lampname, deckname, etc.: s cofraw cafraw
- Change cofraw, cafraw, lampname, deckname, etc.: m lampname =
- Useful link with extra HIRES info:
 - <http://www2.keck.hawaii.edu/inst/hires/startup.html>

Useful directories:

- snaps images are 43" x 57" and have 0.15"/pix for 1x1 binning.
 - 1x1 binning: snap size is:[384, 288] pixels
 - 2x2 binning: snap size is [192, 144] pixels

- 3x3 binning: snap size is [128,96] pixels
- This is for OLD guider images. Needs updated.
- Snaps are transferred nightly to /mir3/snap/fits_todo/YYYYMMDD/.

ThAr Focus Linelist

Output from foc.pro on 9 Jan 2007, for D5, ThAr1, ng3, 4-sec exposure:

Column	Row	Peak	Cts	FWHM	FW@10%	ASYM	Delta	Col
1909	1372	6348	2.22	4.46	-0.15	4.53		
590	1309	519	-Inf	14.00	*****			
2790	1334	1603	2.02	3.98	-0.03	-1.87		
3095	1338	828	2.16	4.32	-0.11	-2.56		
3323	1341	8001	2.29	4.32	-0.13	-2.60		
3705	1347	4007	2.69	5.04	-0.32	-2.99		
478	1260	273	2.75	8.22	-3.23	2.65		
2148	1279	58415	2.58	4.28	0.00	-0.23		
766	1218	17110	2.64	4.92	0.30	0.89		
2271	1235	7144	2.01	3.86	-0.11	0.21		
2968	1244	4744	2.13	4.04	-0.07	-1.66		
2833	1198	375	2.33	6.74	0.93	-0.94		
3043	1200	5020	2.19	4.02	0.00	-3.14		
1490	1139	7276	2.45	4.38	-0.07	0.75		
2010	1145	4382	2.21	4.46	-0.05	0.43		
2713	1153	10029	2.06	3.90	0.06	-1.44		
2485	1109	27365	2.11	3.98	0.03	-1.21		
2284	1067	9564	2.16	4.08	-0.08	0.21		
3031	1076	6247	2.17	4.04	-0.01	-2.20		
1061	1015	25039	2.79	5.04	0.07	0.36		
1852	1024	3250	2.49	4.78	-0.25	0.07		
1380	981	6828	2.58	4.58	0.04	4.12		
1926	987	701	2.33	5.06	1.54	0.33		
3232	1001	15965	2.28	4.30	-0.13	-2.73		
1178	943	780	2.60	6.66	1.32	1.97		
3223	964	14145	2.28	4.24	-0.17	-2.76		
1497	911	1186	2.48	4.76	0.00	1.21		
1940	915	558	2.32	4.28	-0.21	0.18		
2181	918	7569	0.02	0.00	0.00	3.01		
2932	926	61596Sat	2.29	4.32	-0.02	-2.53		
947	871	12924	2.85	5.24	0.20	1.67		
1780	879	11996	2.43	4.48	-0.09	2.36		
2634	888	61561Sat	2.27	4.26	-0.02	-1.62		
1183	1315	719	2.87	5.56	-0.05	0.94		
955	1220	1096	2.56	4.60	0.07	0.38		
3538	1005	959	2.45	4.62	-0.16	-4.54		
2558	1070	3028	2.06	4.00	-0.04	-2.79		